

Phy & Vince

YEE-G02 Samples to reprocess

Date: March 10, 2017

Get well

C=0 C=3

	1	2	3	4	5	6	7	8	9	10	11	12
A	S11 22	S12 22	S13 22	S14 22	S15 22	S16 22	X	X	X	M10 E0 R1	M10 E0 R1	X
B	S11 0	S12 0	S13 0	S14 0	S15 0	S16 0	X	S2 0	X	M10 EU R2	M10 EU R2	X
C	S11 10	S12 10	S13 10	S14 10	S15 10	S16 10	X	S2 10	X	M10 E0 R3	M10 EU R3	X
D	S11 15	S12 15	S13 15	S14 15	S15 15	S16 15	X	S2 15	X	M10 EU R4	M10 E0 R4	X
E	S11 20	S12 20	S13 20	S14 20	S15 20	S16 20	X	S2 20	X	M10 E25 R1	M10 E25 R1	X
F	X	X	X	X	X	X	X	S2 22	X	M10 E25 R2	M10 E25 R2	X
G	S11 24	S12 24	S13 24	S14 24	S15 24	S16 24	X	S2 24	X	M10 E25 R3	M10 E25 R3	X
H	S11 25	S12 25	S13 25	S14 25	S15 25	S16 25	X	S2 25	X	M10 E25 R4	M10 E25 R4	S1 25

Copyright © 2009 Edita Aksamitiene

41

24

Evolution Samples (N= 50; contains all pcr2 samples and all of sample set S2)										
	Sample Set	Replicate	Evolutionary	Competition	a.primer	p1.primer	bc.mean.con	bc.median.con	pcr2	
	S1	NA	25	NA	lon_17	lon_R7	14.1152696	3.38029571	yes	
	S2	NA	0	NA	lon_18	lon_R1	0.08827349	0.06183735	no	
	S2	NA	10	NA	lon_18	lon_R2	1.36412557	0.84897371	no	
	S2	NA	15	NA	lon_18	lon_R3	24.8615259	14.1323601	no	
	S2	NA	20	NA	lon_18	lon_R4	11.4260136	3.65388631	no	
	S2	NA	22	NA	lon_18	lon_R5	20.6938972	20.6820262	yes	
	S2	NA	24	NA	lon_18	lon_R6	22.022152	18.8304087	yes	
	S2	NA	25	NA	lon_18	lon_R7	14.5133095	5.39488153	no	
	S11	NA	0	NA	lon_27	lon_R1	14.1512697	10.6392496	yes	
	S11	NA	10	NA	lon_27	lon_R2	31.4590158	27.8609729	yes	
	S11	NA	15	NA	lon_27	lon_R3	34.1867453	30.8105911	yes	
	S11	NA	20	NA	lon_27	lon_R4	34.5891969	31.6449058	yes	
	S11	NA	22	NA	lon_27	lon_R5	35.1246987	28.4033952	yes	
	S11	NA	24	NA	lon_27	lon_R6	28.7353292	23.6471234	yes	
	S11	NA	25	NA	lon_27	lon_R7	37.8129457	32.8645298	yes	
	S12	NA	0	NA	lon_28	lon_R1	23.3023631	22.1392882	yes	
	S12	NA	10	NA	lon_28	lon_R2	34.8426347	36.4293855	yes	
	S12	NA	15	NA	lon_28	lon_R3	41.5139216	42.1978351	yes	
	S12	NA	20	NA	lon_28	lon_R4	46.1019639	46.3949564	yes	
	S12	NA	22	NA	lon_28	lon_R5	50.1842168	53.161934	yes	
	S12	NA	24	NA	lon_28	lon_R6	23.0165099	21.1650822	yes	
	S12	NA	25	NA	lon_28	lon_R7	44.9864086	44.1823343	yes	
	S13	NA	0	NA	lon_29	lon_R1	29.3383551	29.9844505	yes	
	S13	NA	10	NA	lon_29	lon_R2	42.9601739	36.4727831	yes	
	S13	NA	15	NA	lon_29	lon_R3	36.7913415	31.9589828	yes	
	S13	NA	20	NA	lon_29	lon_R4	44.591743	42.0532527	yes	
	S13	NA	22	NA	lon_29	lon_R5	41.4167466	37.727021	yes	
	S13	NA	24	NA	lon_29	lon_R6	39.6471882	35.4367726	yes	
	S13	NA	25	NA	lon_29	lon_R7	42.6642723	43.0782999	yes	
	S14	NA	0	NA	lon_30	lon_R1	29.9710825	30.4303945	yes	
	S14	NA	10	NA	lon_30	lon_R2	28.7479618	24.261063	yes	
	S14	NA	15	NA	lon_30	lon_R3	34.3190193	26.873492	yes	
	S14	NA	20	NA	lon_30	lon_R4	43.0559351	36.6878952	yes	
	S14	NA	22	NA	lon_30	lon_R5	39.7980915	37.9217795	yes	
	S14	NA	24	NA	lon_30	lon_R6	36.2222689	31.2068756	yes	
	S14	NA	25	NA	lon_30	lon_R7	39.3524257	34.8629139	yes	
	S15	NA	0	NA	lon_31	lon_R1	52.8268995	57.9990636	yes	
	S15	NA	10	NA	lon_31	lon_R2	44.6068871	44.750351	yes	
	S15	NA	15	NA	lon_31	lon_R3	21.0522351	16.9405569	yes	
	S15	NA	20	NA	lon_31	lon_R4	48.167732	52.6620179	yes	
	S15	NA	22	NA	lon_31	lon_R5	47.9421064	47.3551564	yes	
	S15	NA	24	NA	lon_31	lon_R6	48.1970821	50.4575971	yes	
	S15	NA	25	NA	lon_31	lon_R7	35.3397788	31.7390849	yes	
	S16	NA	0	NA	lon_32	lon_R1	26.9911755	19.7845528	yes	
	S16	NA	10	NA	lon_32	lon_R2	58.755851	63.2331152	yes	
	S16	NA	15	NA	lon_32	lon_R3	21.8175747	17.2786558	yes	
	S16	NA	20	NA	lon_32	lon_R4	43.9097998	39.8296687	yes	
	S16	NA	22	NA	lon_32	lon_R5	20.6203919	15.0336288	yes	
	S16	NA	24	NA	lon_32	lon_R6	56.5033875	49.824418	yes	
	S16	NA	25	NA	lon_32	lon_R7	41.9039173	33.3038751	yes	

Competition Samples (N= 16; contains all of mix M10)

1	M10	R1	0	0	0 lon_33	lon_R11	19.905329	19.158917	no
2	M10	R2	0	0 lon_34	lon_R11	8.818273	7.671685	no	
3	M10	R3	0	0 lon_35	lon_R11	2.357958	2.128046	no	
4	M10	R4	0	0 lon_36	lon_R11	2.597363	2.184959	no	
5	M10	R1	25	0 lon_37	lon_R11	18.183827	11.227402	no	
6	M10	R2	25	0 lon_38	lon_R11	50.161297	48.882898	no	
7	M10	R3	25	0 lon_39	lon_R11	23.629128	11.153746	no	
8	M10	R4	25	0 lon_40	lon_R11	33.003047	27.801369	no	
9	M10	R1	0	3 lon_25	lon_R12	35.046062	29.249595	no	
10	M10	R2	0	3 lon_26	lon_R12	52.713971	56.093657	no	
11	M10	R3	0	3 lon_27	lon_R12	30.290724	28.500221	no	
12	M10	R4	0	3 lon_28	lon_R12	10.544308	10.073059	no	
13	M10	R1	25	3 lon_29	lon_R12	41.028604	36.21226	no	
14	M10	R2	25	3 lon_30	lon_R12	3.550894	1.673805	no	
15	M10	R3	25	3 lon_31	lon_R12	32.16729	26.177167	no	
16	M10	R4	25	3 lon_32	lon_R12	4.870915	2.431687	no	

Layout for DNA extraction

Plate T=0	
A	T0.M1
B	T0.M2
C	T0.M1
D	T0.M1
E	T25.M1
F	T25.M2
G	T25.M1
H	T25.M1

Plate T=3	
A	T0.M2
B	T0.M2
C	T0.M2
D	T0.M1
E	T25.M1
F	T25.M2
G	T25.M1
H	T25.M1

Layout for Primer addition

Plate T=0	
Ion_33	1
Ion_34	2
Ion_35	3
Ion_36	4
Ion_37	5
Ion_38	6
Ion_39	7
Ion_40	8

Plate T=3	
Ion R8	97
Ion R9	98
Ion R10	99
Ion R11	100
Ion R12	101
Ion R13	102

Layout for Ion R

Plate T=0	
Ion R1	1
Ion R2	2
Ion R3	3
Ion R4	4
Ion R5	5
Ion R6	6
Ion R7	7
Ion R8	8
Ion R9	9
Ion R10	10
Ion R11	11
Ion R12	12
Ion R13	13

Plate T=3 (cor Ion R)	
Ion R8	145
Ion R9	146
Ion R10	147
Ion R11	148
Ion R12	149
Ion R13	150

Plate T=3	
Ion_17	97
Ion_18	98
Ion_19	99
Ion_20	100
Ion_21	101
Ion_22	102
Ion_23	103
Ion_24	104

	2 Do PCR using A and P1 primers for DNA plate 1 and 2	PCR plate label: FIT PCR 4-1, 4-2
a	DNA (lyticase+beads)	ul
	primer A+P1	0.25+0.25
	Taq plus mix (Lamda Cat#D124R)	12.5
	water	7.5
	total volume	25
b	program: UPTAG	
	94C	5min
	94C	10sec
	50C	35 cycles
	72C	10sec
	72C	10sec
c	run 2% gel to check PCR result	5min
		clean single band of 131bp

	Ion_17	Ion_18	Ion_19	Ion_20	Ion_21	Ion_22	Ion_23	Ion_24	Ion_25	x	x	x
Ion_R1	x	x	x	x	x	x	S1_25	M10_E0_C0	M10_E0_C3	x	x	x
Ion_R2	S11_0	S12_0	S13_0	S14_0	S15_0	S16_0	S2_0	M10_E0_C0	M10_E0_C3	x	x	x
Ion_R3	S11_10	S12_10	S13_10	S14_10	S15_10	S16_10	S2_10	M10_E0_C0	M10_E0_C3	x	x	x
Ion_R4	S11_15	S12_15	S13_15	S14_15	S15_15	S16_15	S2_15	M10_E0_C0	M10_E0_C3	x	x	x
Ion_R5	S11_20	S12_20	S13_20	S14_20	S15_20	S16_20	S2_20	M10_E25_C0	M10_E25_C3	x	x	x
Ion_R6	S11_22	S12_22	S13_22	S14_22	S15_22	S16_22	S2_22	M10_E25_C0	M10_E25_C3	x	x	x
Ion_R7	S11_24	S12_24	S13_24	S14_24	S15_24	S16_24	S2_24	M10_E25_C0	M10_E25_C3	x	x	x
Ion_R8	S11_25	S12_25	S13_25	S14_25	S15_25	S16_25	S2_25	M10_E25_C0	M10_E25_C3	x	x	x

l z 3 4 5 6 7 8 9

Protocol: Lyticase + beads DNA extraction

Date Edited: April 17, 2017

1. Prepare the samples
 - a. Add 60ul of yeast sample in each appropriate well of a 96-well pcr microplate
 - b. Centrifuge at 3000 rcf for 5 minutes to pellet
 - c. Pull the supernatant one well at a time, don't need to get all of it.
 - d. Store covered in freezer (-20 C) if the rest of the protocol will not be run immediately
2. Prepare the lyticase and beads mixture
 - a. Fill a single microcentrifuge tube about 1/3 of the way with lyticase solution (100 units / ml).
 - b. Add glass beads (0.5 mm diameter) to the tube for a total volume of approximately 2/3 of the microcentrifuge tube's capacity.
 - c. Refill the tube with beads or lyticase solution as needed throughout the protocol. Keep the lyticase stock tube on ice and the lyticase + beads mixture on cooled plate holder when not in use.
3. Add the lyticase and beads mixture to each well
 - a. Cut the tip off of a large pipette tip (p1000 size) at a 45 degree angle to allow both lyticase and beads to be drawn into the tip.
 - b. Transfer 100ul of lyticase + beads to each well (do not touch the well or sample, change pipette tips if you do) – aim for ~75% lyticase and ~25% beads per well.
4. Resuspend / bead-beat the samples before PCR
 - a. Cover the plate with a microseal 'B' seal cover (brown backing) completely, ensuring that each individual well is completely sealed so that there is no contamination between wells.
 - b. Vortex the plate at the highest setting
 - i. Push down on each pair of columns for approximately 2 seconds of vortexing to ensure that all wells are mixed evenly and completely
 - ii. Repeat as needed until all wells are resuspended
5. Incubate the PCR plate in the thermocycler
 - a. Load the plate into the pcr machine ensuring that it is firmly in place
 - b. Navigate to the program named "37C_30" on the thermocycler and double check that the program specifications match those outlined in appendix 1.
 - c. close the thermocycler lid and press run.
 - i. Do not use the rubberish compression pad that is stored near the vortexer.
 - d. Pause the PCR program after 15 minutes and vortex / bead-beat all samples again as specified in the "resuspend" step above.
 - i. Double check before vortexing that all wells are still fully sealed by the pcr cover to avoid contamination across wells.
 - e. Place the plate back in the thermocycler and resume the PCR program. Allow the program to run to completion.
 - i. *If tubes are allowed to transition fully from the highest temperature at the end of the program to the 4 degrees C holding temperature they will be highly deformed. Pulling the tubes out as soon as the program has finished running will prevent this deformation.
6. Spin down the plate (Hermle Z 300)
 - a. Spin at 3000 rcf for 5 mminutes
 - b. Remove supernatant (50ul per well) and add to a new pcr-plate ("DNA storage plate")
 - i. Will only need 5ul for each round of pcr, so there should be a lot of sample in reserve for future use
 - ii. OK to use multichannel pipette for this step
 - c. If you are stopping here, the plate can be stored at -20 C in the freezer until you are ready to do the PCR protocol.

Appendix 1: Thermocycler program specifications

1. Parameters
 - a. Volume = 100 μ l
 - b. Lid = 105 C
2. Protocol
 - a. 37 C for 30 minutes
 - b. 95 C for 10 minutes
 - c. 4 C for infinite hold

	Ion_17	Ion_18	Ion_19	Ion_20	Ion_21	Ion_22	Ion_23	Ion_24	Ion_25	x	x	x	x
Ion_R1	x	x	x	x	x	x	S15_0	S15_25	M10_E0_C0	M10_E0_C3	x	x	x
Ion_R2	S12_0	S13_0	S14_0	S15_0	S16_0	S2_0	M10_E0_C0	M10_E0_C3	x	x	x	x	x
Ion_R3	S12_10	S13_10	S14_10	S15_10	S16_10	S2_10	M10_E0_C0	M10_E0_C3	x	x	x	x	x
Ion_R4	S12_15	S13_15	S14_15	S15_15	S16_15	S2_15	M10_E0_C0	M10_E0_C3	x	x	x	x	x
Ion_R5	S12_20	S13_20	S14_20	S15_20	S16_20	S2_20	M10_E25_C0	M10_E25_C3	x	x	x	x	x
Ion_R6	S12_22	S13_22	S14_22	S15_22	S16_22	S2_22	M10_E25_C0	M10_E25_C3	x	x	x	x	x
Ion_R7	S12_24	S13_24	S14_24	S15_24	S16_24	S2_24	M10_E25_C0	M10_E25_C3	x	x	x	x	x
Ion_R8	S12_25	S13_25	S14_25	S15_25	S16_25	S2_25	M10_E25_C0	M10_E25_C3	x	x	x	x	x

	Ion_17	Ion_18	Ion_19	Ion_20	Ion_21	Ion_22	Ion_23	Ion_24	Ion_25	x	x	x	x
Ion_R1	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Ion_R2	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Ion_R3	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Ion_R4	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Ion_R5	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Ion_R6	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Ion_R7	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Ion_R8	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

	Group 1	Group 2
Vol to add	Total volume	Total volume

→ S12 - 22	7.86	13.99	V	154.39 . ml	(175)
S12 - 24	9.16	12.01	V	154.39 . ml	(175)
S12 - 25	8.02	13.72	V	154.39 . ml	(175)
S13 - 0	5.70	19.30	V	154.39 . ml	(175)
S13 - 10	9.58	11.48	V	266.60 ml (1,363)	(500 ml)
S14 - 15	7.76	14.16	V	500 ml	(500 ml)

Group 1 : 20 samples
Group 2 : 46 samples

Group 2
Group 1

Group 2
Group 1

Group 2
Group 1

Group 2
Group 1

↳ Increase Sample volume relative to water volume in PCR protocol to sample DNA yields

Group 1 : 20 samples
Group 2 : 46 samples

Sample ID	N Wells	N Barcodes	A Primer	P1 Primer	PCR (ng/uL)	PCR (ng/uL) [RUN 2]	Vol to add	Vol to add rerun	Total Volume	Total volume for pool (rerun)
S11_0	11	22	Ion_17	Ion_R2	25.2		4.37		1	
S11_10	11	22	Ion_17	Ion_R3	19.4		5.67		1	
S11_15	11	22	Ion_17	Ion_R4	17.4		6.32		1	
S11_20	11	22	Ion_17	Ion_R5	13.2		8.33		1	
S11_22	11	22	Ion_17	Ion_R6	20.8		5.29		1	
S11_24	11	22	Ion_17	Ion_R7	17.9		6.15		1	
S11_25	11	22	Ion_17	Ion_R8	17.3		6.36		1	
S12_0	11	22	Ion_18	Ion_R2	19.6		5.61		1	
S12_10	11	22	Ion_18	Ion_R3	20.4		5.39		1	
S12_15	11	22	Ion_18	Ion_R4	19.2		5.73		1	
S12_20	11	22	Ion_18	Ion_R5	15		7.33		1	
S12_22	11	22	Ion_18	Ion_R6	3.32	7.86	33.13	13.99	1	
S12_24	11	22	Ion_18	Ion_R7	3.26	9.16	33.74	12.01	1	
S12_25	11	22	Ion_18	Ion_R8	4.28	8.02	25.70	13.72	1	
S13_0	11	22	Ion_19	Ion_R2	3	5.7	36.67	19.30	1	
S13_10	11	22	Ion_19	Ion_R3	2.4	9.58	45.83	11.48	1	
S13_15	11	22	Ion_19	Ion_R4	21.6		5.09		1	
S13_20	11	22	Ion_19	Ion_R5	20.8		5.29		1	
S13_22	11	22	Ion_19	Ion_R6	14.6		7.53		1	
S13_24	11	22	Ion_19	Ion_R7	18.1		6.08		1	
S13_25	11	22	Ion_19	Ion_R8	15.1		7.28		1	272.90
S14_0	11	22	Ion_20	Ion_R2	17.2		6.40		2	
S14_10	11	22	Ion_20	Ion_R3	21		5.24		2	
S14_15	11	22	Ion_20	Ion_R4	2.68	7.76	41.04	14.18	2	
S14_20	11	22	Ion_20	Ion_R5	7.12		15.45		2	
S14_22	11	22	Ion_20	Ion_R6	16		6.88		2	
S14_24	11	22	Ion_20	Ion_R7	13		8.46		2	
S14_25	11	22	Ion_20	Ion_R8	15.6		7.05		2	
S15_0	11	22	Ion_21	Ion_R2	17.3		6.36		2	
S15_10	11	22	Ion_21	Ion_R3	13.1		8.40		2	
S15_15	11	22	Ion_21	Ion_R4	6.76		16.27		2	
S15_20	11	22	Ion_21	Ion_R5	16.4		6.71		2	
S15_22	11	22	Ion_21	Ion_R6	16.8		6.55		2	
S15_24	11	22	Ion_21	Ion_R7	12.2		9.02		2	
S15_25	11	22	Ion_21	Ion_R8	15.3		7.19		2	
S16_0	11	22	Ion_22	Ion_R2	16.1		6.83		2	
S16_10	11	22	Ion_22	Ion_R3	15.2		7.24		2	
S16_15	11	22	Ion_22	Ion_R4	11.7		9.40		2	
S16_20	11	22	Ion_22	Ion_R5	17.4		6.32		2	
S16_22	11	22	Ion_22	Ion_R6	17.6		6.25		2	
S16_24	11	22	Ion_22	Ion_R7	12.1		9.09		2	
S16_25	11	22	Ion_22	Ion_R8	10.6		10.38		2	
S1_25	6	12	Ion_23	Ion_R1	15.8		3.80		2	
S2_0	5	10	Ion_23	Ion_R2	16.3		3.07		2	
S2_10	5	10	Ion_23	Ion_R3	2.9		17.24		2	
S2_15	5	10	Ion_23	Ion_R4	17.8		2.81		2	
S2_20	5	10	Ion_23	Ion_R5	12		4.17		2	
S2_22	5	10	Ion_23	Ion_R6	13.6		3.68		2	
S2_24	5	10	Ion_23	Ion_R7	3.18		15.72		2	
S2_25	5	10	Ion_23	Ion_R8	15.9		3.14		2	
M10_E0_CO_R1	1	12 + ref @50%	Ion_24	Ion_R1	14.1		0.71		2	
M10_E0_CO_R2	1	12 + ref @50%	Ion_24	Ion_R2	7.32		1.37		2	
M10_E0_CO_R3	1	12 + ref @50%	Ion_24	Ion_R3	11.2		0.89		2	
M10_E0_CO_R4	1	12 + ref @50%	Ion_24	Ion_R4	13.8		0.72		2	
M10_E25_CO_R1	1	12 + ref @50%	Ion_24	Ion_R5	11.2		0.89		2	
M10_E25_CO_R2	1	12 + ref @50%	Ion_24	Ion_R6	3.74		2.67		2	
M10_E25_CO_R3	1	12 + ref @50%	Ion_24	Ion_R7	14.9		0.67		2	
M10_E25_CO_R4	1	12 + ref @50%	Ion_24	Ion_R8	13.5		0.74		2	
M10_E0_C3_R1	1	12 + ref @50%	Ion_25	Ion_R1	4.2		2.38		2	
M10_E0_C3_R2	1	12 + ref @50%	Ion_25	Ion_R2	5.8		1.72		2	
M10_E0_C3_R3	1	12 + ref @50%	Ion_25	Ion_R3	16.7		0.60		2	
M10_E0_C3_R4	1	12 + ref @50%	Ion_25	Ion_R4	12.3		0.81		2	
M10_E25_C3_R1	1	12 + ref @50%	Ion_25	Ion_R5	3.84		2.60		2	
M10_E25_C3_R2	1	12 + ref @50%	Ion_25	Ion_R6	16.9		0.59		2	
M10_E25_C3_R3	1	12 + ref @50%	Ion_25	Ion_R7	14.8		0.68		2	
M10_E25_C3_R4	1	12 + ref @50%	Ion_25	Ion_R8	7.84		1.28		2	273.08
							552.38		252.61	

GROUP	sample/grou	concentration	volume to mix	vol to mix / 2
G1	220	7.78	28.277635	14.14
G2	299	10.2	29.3137255	14.66

Water to add

519 H2O: 57.5913605 28.80







add up 1ng / 519 129.75 64.875 36.08

F2-2 131bp PCR products 4 ng/ul

50nM

≤ 20,000 counts = 
 Contamination ≥ 5% = 

	Ion_17	Ion_18	Ion_19	Ion_20	Ion_21	Ion_22	Ion_23	Ion_24	Ion_25	Ion_26	x	x
Ion_R1	x	x	x	x	x	x	S1	25	M10_E0_C0	M10_E0_C3	x	x
Ion_R2	S12_0	S13_0	S14_0	S15_0	S16_0	S2_0	S2_10	S2_15	M10_E0_C0	M10_E0_C3	x	x
Ion_R3	S12_10	S13_10	S14_10	S15_10	S16_10	S2_10	S2_15	S2_20	M10_E0_C0	M10_E0_C3	x	x
Ion_R4	S12_15	S13_15	S14_15	S15_15	S16_15	S2_15	S2_20	S2_22	M10_E25_C0	M10_E25_C1	x	x
Ion_R5	S11_20	S12_20	S13_20	S14_20	S15_20	S16_20	S16_22	S16_24	M10_E25_C0	M10_E25_C1	x	x
Ion_R6	S12_22	S13_22	S14_22	S15_22	S16_22	S2_22	S2_24	S2_25	M10_E25_C0	M10_E25_C1	x	x
Ion_R7	S12_24	S13_24	S14_24	S15_24	S16_24	S2_24	S2_25	S2_25	M10_E25_C0	M10_E25_C1	x	x
Ion_R8	S12_25	S13_25	S14_25	S15_25	S16_25	S2_25	S2_25	S2_25	M10_E25_C0	M10_E25_C1	x	x

	Ion_17	Ion_18	Ion_19	Ion_20	Ion_21	Ion_22	Ion_23	Ion_24	Ion_25	Ion_26	x	x
Ion_R1	4.37	5.61	36.67	6.40	6.36	6.83	3.80	0.71	2.38			
Ion_R2	5.67	5.39	45.83	5.24	8.40	7.24	3.07	1.37	1.72			
Ion_R3	6.32	5.73	5.09	41.04	16.27	9.40	17.24	0.89	0.60			
Ion_R4	8.33	7.33	5.29	15.45	6.71	6.32	2.81	0.72	0.81			
Ion_R5	5.29	33.13	7.53	6.88	6.55	6.25	4.17	0.89	2.60			
Ion_R6	6.15	33.74	6.08	8.46	9.02	9.09	15.72	0.67	0.59			
Ion_R7	6.36	25.70	7.28	7.05	7.19	10.38	3.14	0.74	1.28			

Group 1	Group 2
---------	---------

